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EN
MANUFACTURER'S INSTRUCTIONS AND INFORMATION - INSTRUCTIONS
FOR USE OF THE PERSONAL PROTECTIVE EQUIPMENT

easyCARE nitrile Nitrile examination gloves, powder-free

REF: RNBRRPPP01

RR - size, PPP - number of pieces per box

Device features

- raw material: 100% nitrile (acrylonitrile butadiene)
- colour: blue
- cuff: beaded
- powder: powder-free, the content of powder residue below 2mg per glove
- AQL: 1,0
- internal surface: chlorinated and polymerized
- outer surface: polymerized, micro-textured throughout the glove with additional texture at the fingertips
- allowed for contact with food: in accordance with Regulation (EC) 1935/2004 and Commission Regulation (EC) 10/2011
- intermediate packaging: 100 pcs or 200 pcs with a colour-coded packaging for easy size recognition
- size: XS-XL
- single-use
- non-sterile
- shelf life: 5 years

Intended use

The device is intended for patient diagnostics, providing a barrier to the transmission of microorganisms between the operator and the patient, minimizing the risk of cross-contamination.

Product classification

- Class I Medical Device according to MDR 2017/745 (EU)
- Personal Protective Equipment Cat.III according to Regulation (EU) 2016/425 EU Type Examination and Ongoing Conformity by Notify Body: SATRA Technology Europe Ltd, Bracetown Business Park Clonee, D15 YN2P, Ireland (Notified Body: 2777)

Compliance with standards:

EN ISO 21420:2020, EN ISO 21420:2020/A1:2024, EN ISO 374-1:2016, EN ISO 374-1:2016/A1:2018, EN ISO 374-2:2019, EN 16523-1:2015+A1:2018, EN ISO 374-4:2019, EN ISO 374-5:2016, EN 455-1:2020+A2:2024, EN 455-2:2024, EN 455-3:2023, EN 455-4:2009, ASTM F1671, ASTM D6978, ASTM D6319

Contraindications

The ingredients used to make the gloves may cause allergic reactions in some people. If these occur, discontinue use and contact a medical specialist.

Research conducted:

Chemical penetration (EN ISO 374-1:2016+A1:2018/Type B)	Level	Mean Degradation (EN ISO 374-4:2019)
K 40% Sodium hydroxide	6	-33,5%
P 30% Hydrogen peroxide	4	22,1%
T 37% Formaldehyde	6	-0,3%

Degradation levels indicate the change in puncture resistance after exposure to the challenge chemical.

10%	Ethanol	6
50%	Benzalkonium chloride	6
70%	Citric acid	6
10%	Hydrochloric acid	6
10-13%	Sodium hypochlorite	6
35%	Sodium nitrate	6
-	Ethidium bromide (saturated)	6
4%	Glutaraldehyde	6

Levels of resistance to penetration:

Level 1: >10 min; Level 2: >30 min; Level 3: >60 min;
Level 4: >120 min; Level 5: >240 min; Level 6: >480 min.

Testing according to ASTM D 6978 Test Chemotherapy Drugs and Concentration	Minimum breakthrough time Detection time (min)
Carmustine (BCNU) 3.3 mg/ml (3,300ppm)	24.1
Cisplatin 1.0 mg/ml (1,000 ppm)	240
Cyclophosphamide (Cytoxan) 20 mg/ml (20,000ppm)	240
Doxorubicin Hydrochloride 2.0 mg/ml (2,000ppm)	240
Etoposide (Toptsar) 20.0 mg/ml (20,000ppm)	240
Fluorouracil 50.0 mg/ml (50,000ppm)	240
Mechlorethamine HCl 1.0 mg/ml (1,000ppm)	240
Methotrexate 25 mg/ml (25,000ppm)	240
Mitomycin C 0.5 mg/ml (500ppm)	240
Paclitaxel (Taxol) 6.0 mg/ml (6,000ppm)	240
Thiotepa 10.0 mg/ml (10,000ppm)	64.3
Vincristine Sulfate 1.0 mg/ml (1,000ppm)	240

Resistance against Bacteria and Fungi - PASS, according to EN ISO 374-5:2016

Resistance against Virus - PASS, according to EN ISO 374-5:2016

User information:

1. Application:

These gloves are intended for single use only (up to 60 min.), on one individual during a single procedure, and should be replaced with new gloves after use.

2. Duration of use:

The information on the duration of use does not reflect the actual duration of protection in the workplace and the differentiation between mixtures and pure chemicals.

3. Chemical resistance:

The chemical and penetration resistance of the glove has been assessed under laboratory conditions from samples taken from the palm only, and relates only to the chemical tested (this may be different if the chemical is used in a mixture).


4. Workplace protection:

- It is recommended to check that the gloves are suitable for the intended use because the conditions at the workplace may differ from the type test depending on temperature, abrasion and degradation.

- During use, protective gloves may provide less resistance to the dangerous chemical due to changes in physical properties. Movements, snagging, rubbing, degradation caused by the chemical contact etc. may reduce the actual use time significantly. For corrosive chemicals, degradation can be the most important factor to consider in selecting chemical-resistant gloves.

- Before use, check the gloves for defects or damage and that the packaging is not damaged (torn).

Putting on and removing gloves:

Sizes:

X-Small	5-6
Small	6-7
Medium	7-8
Large	8-9
X-Large	9-10
XX-Large	10-11

Additional information:

- Mechanically packed in layers in a package.
- For more information (including Declaration of Conformity) visit our website: www.zarys.com
- easyCARE gloves protect against hazardous substances and mixtures as well as harmful biological agents.

Special storage or handling conditions

- The product has a shelf life of 5 years from the date of manufacture.
- Store in a cool, dry, well-ventilated and clean place at temperature: 5°C–35°C away from direct sunlight. Keep away from sources of ozone or fire. Store the gloves in the original packaging. Improper storage reduces the protective properties of the gloves.

Disposal

Used devices may pose a potential biohazard. Devices should be disposed of in accordance with local medical waste regulations.

The easyCARE gloves are classified as a Class I medical device in accordance with Regulation 2017/745 (Annex VIII) and a Category III PPE in accordance with Regulation (EU) 2016/425 (Annex I), however, this document does not constitute instructions for use of the medical device. As a Class I medical device, the device can be safely used without instructions for use and based on the provisions of Annex I p. 23.1(d) of Regulation (EU) 2017/745, instructions for use for the medical device are not required. At the same time, these instructions in no way constitute instructions for use for a medical device and must not be used for that purpose.

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