


# Safeflow Extension Set

|           |                           |           |
|-----------|---------------------------|-----------|
| <b>en</b> | Instructions for use      | <b>3</b>  |
| <b>de</b> | Gebrauchsanweisung        | <b>4</b>  |
| <b>bg</b> | Инструкция за употреба    | <b>5</b>  |
| <b>cs</b> | Návod k použití           | <b>7</b>  |
| <b>da</b> | Brugsanvisning            | <b>8</b>  |
| <b>el</b> | Οδηγίες χρήσης            | <b>9</b>  |
| <b>es</b> | Instrucciones de uso      | <b>10</b> |
| <b>et</b> | Kasutusjuhised            | <b>12</b> |
| <b>fi</b> | Käyttöohje                | <b>13</b> |
| <b>fr</b> | Mode d'emploi             | <b>14</b> |
| <b>hr</b> | Upute za upotrebu         | <b>15</b> |
| <b>hu</b> | Használati utasítás       | <b>16</b> |
| <b>id</b> | Petunjuk penggunaan       | <b>18</b> |
| <b>it</b> | Istruzioni per l'uso      | <b>19</b> |
| <b>lt</b> | Naudojimo instrukcija     | <b>21</b> |
| <b>lv</b> | Lietošanas instrukcija    | <b>22</b> |
| <b>nl</b> | Gebruikersinformatie      | <b>23</b> |
| <b>no</b> | Bruksanvisning            | <b>24</b> |
| <b>pl</b> | Instrukcja użytkowania    | <b>25</b> |
| <b>pt</b> | Instruções de utilização  | <b>27</b> |
| <b>ro</b> | Instrucțiuni de utilizare | <b>28</b> |
| <b>ru</b> | Инструкция по применению  | <b>30</b> |
| <b>sk</b> | Návod na použitie         | <b>31</b> |
| <b>sl</b> | Navodila za uporabo       | <b>32</b> |
| <b>sr</b> | Uputstvo za upotrebu      | <b>34</b> |
| <b>sv</b> | Bruksanvisning            | <b>35</b> |
| <b>th</b> | คำแนะนำในการใช้งาน        | <b>36</b> |
| <b>vi</b> | hướng dẫn sử dụng         | <b>37</b> |
| <b>zh</b> | 使用说明                      | <b>38</b> |

**B | BRAUN**

## **en** Instructions for use

### **Device description**

Safeflow is a luer activated valve that is developed as needlefree injection port in IV applications. Safeflow Extension Set is used to extend infusion lines by attaching them to a luer lock connection.

Safeflow Extension Sets can be used with devices which are in conformance with Luer connectors (e.g. ISO 80369-7). For gravity or pressure use up to 2 bar.

### **Sterility**

Sterilized using ethylene oxide.

Please refer to the label of primary packaging.

### **Intended purpose**

Smallbore extension set for infusion lines.

### **Indication**

**Safeflow Extension Set without back-check valve:** Infusion and transfusion therapy as clinically indicated and according to the SPC of the drugs / solutions.

**Safeflow Extension Sets with back-check valve:** Infusion therapy as clinically indicated and according to the SPC of the drugs / solutions.

### **Patient population**

**Safeflow Extension Set without back-check valve:** Safeflow Extension Sets can be used for all patients for whom infusion or transfusion therapy is prescribed.

**Safeflow Extension Sets with back-check valve:** Safeflow Extension Sets can be used for all patients for whom infusion therapy is prescribed.

### **Intended User**

Safeflow Extension Sets should only be used by authorized persons e.g. healthcare professionals, licensed and/or certified physicians, nurses, pharmacists, caregivers (according to local regulations) who have been adequately trained within the educational background in this technique. Upon medical assessment and adequate instruction, patients and / or caregivers may also be allowed to take over defined handling steps according to national guidelines. Healthcare workers are responsible for ensuring that patients and / or caregivers are instructed in the correct use of the product and that the instruction of caregiver and / or patient is documented by the healthcare worker.

### **Contraindications**

**Safeflow Extension Set with back-check valve** cannot be used for transfusion therapy as there might be the risk that blood clots may form in the body of the back-check valve.

### **Residual Risks/Side effects**

General risks occurring during infusion or transfusion therapy include air embolism, local and systemic infection, particle contamination, medication error including over and under supply. Kinking of the line may result in decrease or stop of infusion flow. Leakage of drugs / fluids can be an additional hazard to patients, visitors or staff. Transfusion-specific side effects are possible and needs to be monitored.

### **Precautions**

- Visually inspect the sterile barrier system for breaches and integrity prior to use.
- Please also consider the corresponding SPCs of the drugs/solutions regarding possible incompatibilities of drugs/solutions.
- Prime product by attaching a priming device before usage
- Flush product per national and / or institutional protocol.
- Mind the manufacturer's instructions for use of complementarily used products (e.g. IV administration Set or infusion pump).

### **Warning**

- Do not use if packaging or product is damaged or contaminated or if protective caps are loose or missing. Use of damaged or contaminated device might lead to injury, illness or death of the patient.
- Do not re-use. The product is a single-use, sterile and nonpyrogenic medical device. Re-use of single-use devices creates a potential risk of patient or user. It may lead to contamination and / or impairment of functional capability. Contamination and / or limited functionality of the device may lead to injury, illness or death of the patient.
- Monitor all connections for tightness over the whole application
- Luer slip syringes should be hand-tightened into the valve using a 90° turn. Do not leave luer slip syringes unattended.
- Syringes or connectors which are not in compliance with ISO 80369-7 / EN ISO 80369-7 can damage the swabable valve. Syringes and male luer connectors have a large variety of configurations and can vary significantly in design and dimensions.
- Do not use, in any situation, needles in the needle-free injection port.
- Do not re-attach the protective cap of the patient connector.
- Close clamp before disconnecting the syringe. Turn the syringe counter-clockwise to remove.
- In case of continuous infusion of rapid flow rates, e.g. in case of volume replacement, the presence of needleless connectors may reduce flow rates, in particular when used in combination with large-bore catheters. Negative clinical outcomes might result when therapies with rapid flow rates are impeded.

### **Operating Instructions**

Mind institutional protocols and national guidelines for hygienic procedures such as disinfection and flushing during entire process.

### **Duration of use**

Safeflow Extension Sets should be changed acc. to national standards (e.g. CDC) and / or institutional protocols. Mind that duration of use is dependent on the intended therapy according to the SPC of the solution or drug.

### **Disposal**

Disposal acc. to local guidelines and / or hospital protocols.

### **Storage and handling conditions**

Storage conditions:

Keep away from sunlight.

Keep dry.

### **Notice to the user**

If, during the use of this product or as a result of its use, a serious incident has occurred, please report it to the manufacturer and / or its authorized representative and to your national authority. If further Instructions for Use are required, they can be requested from the manufacturer or can be obtained from the B. Braun homepage: <https://eifu.bbraun.com/>.

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B. Braun Melsungen AG,
B. Braun (Thailand) Ltd.,
B. Braun Medical AG