
Transofix®


en	Instructions for use	4
de	Gebrauchsanweisung	5
bg	Инструкция за употреба	6
cs	Návod k použití	7
da	Brugsanvisning	8
el	Οδηγίες χρήσης	9
es	Instrucciones de uso	10
et	Kasutusjuhised	11
fi	Käyttöohjeet	12
fr	Mode d'emploi	13
hr	Upute za uporabu	14
hu	Használati utasítás	15
id	Petunjuk penggunaan	16
it	Istruzioni per l'uso	17
lt	Naudojimo instrukcija	18
lv	Lietošanas instrukcija	19
nl	Gebruikersinformatie	20
no	Bruksanvisning	21
pl	Instrukcja użytkowania	22
pt	Instruções de utilização	23
ro	Instrucțiuni de utilizare	24
ru	Инструкция по применению	25
sk	Návod na používanie	26
sl	Navodila za uporabo	27
sr	Uputstvo za upotrebu	28
sv	Bruksanvisning	29
zh	使用说明	30

Transofix®

en	Transfer set for sterile fluids
de	Transferset für sterile Flüssigkeiten
bg	Комплект за прехвърляне на стерилни течности
cs	Souprava pro přenos sterilních tekutin
da	Overførselssæt til sterile væsker
el	Σετ μεταφοράς για αποστειρωμένα υγρά
es	Kit para transferencia de líquidos estériles
et	Steriilsete vedelike kandekomplekt
fi	Siirtosarja steriileille nesteille
fr	Set de transfert pour liquides stériles
hr	Set za prijenos sterilnih tekućina
hu	Átviteli készlet steril folyadékokhoz
id	Set transfer untuk cairan steril
it	Set per il trasferimento sterile di liquidi
lt	Sterilių skysčių perpylimo rinkinys
lv	Sterilu šķidrumu pārļiešanas komplekts
nl	Overdrachtset voor steriele vloeistoffen
no	Overføringssett for sterile væsker
pl	Zestaw do przenoszenia jałowych płynów
pt	Sistema de transferencia de Fluidos Estéreis
ro	Set de transfer pentru lichide sterile
ru	Двухсторонняя канюля для смешивания стерильных инфузионных растворов
sk	Prevodová hadička pre sterilné tekutiny
sl	Set za pretok sterilnih tekočin
sr	Komplet za transfer sterilnih tečnosti
sv	Överföringskanyl för sterila lösningar
zh	无菌液体转移接头

en Instructions for use

Device description

Transofix® has an interface to all types of infusion containers with a port system and medication vials.

Transofix® is a transfer set for sterile fluids.

Materials used:

Protective Caps: PE

Piercing Spike: SAN/ABS

Do not re-use.

Sterility

Sterilized using irradiation.

Do not use if package is damaged.

Visually inspect the sterile barrier system for breaches and integrity prior to use.

Intended purpose

Transfer set for sterile fluids.

Indication

- Liquid transfer between different containers.
- Reconstitution of lyophilized drugs according to the Summary of product characteristics (SPC) of the drugs.
- Addition of drugs.
- Dilution of solid and liquid drugs according to the Summary of product characteristics of the drugs.
- Transfer of transfusions.

Patient population

Transofix® is not directly used on patients as it is only used for drug admixture.

Intended User

The products stated above are intended to be used by the following users:

- pharmacy health care professionals
- nursing professionals
- physicians
- caretakers (Homecare group)

Contra-indications

No contraindications.

Side effects

No side effects related to the device are known.

The following complications may occur during the procedure of liquid transfer (drugs, plasma, blood, parenteral nutrition), reconstitution or dilution:

- Microbiological or particulate contamination as a result of poor aseptic technique up to infection
- Medication errors

Precautions

Always refer to the Summary of product characteristics of the drugs/solutions for admixture procedure.

Warning

- Do not use if packaging or product is damaged or contaminated. Use of damaged or contaminated device might lead to injury, illness or death of the patient.
- If product is used with cytotoxic or antineoplastic drugs, please mind safe handling instructions and guidelines/institutional protocols.
- Mind the Summary of product characteristics of the drugs for material compatibility.
- Re-use of single-use devices creates a potential risk of patient or user. It may lead to contamination and/or impairment of functional capability. Contamination and/or limited functionality of the device may lead to injury, illness or death of the patient.
- Do not re-sterilize.

Operating Instructions

1. Remove first protective cap, place the Transofix® on the upright standing container/vial and puncture the port using the large grip flanges, avoiding rotation.
2. Keep the punctured container/vial upright
3. Remove second protective cap, hold the second container/vial in an upside-down position over the free spike and puncture its port

using one strong movement to avoid leakage, avoiding rotation.

4. Transfer the drug/fluid into the chosen container/vial.
5. Finally, remove the empty container/vial together with the Transofix®.

Duration of use

Transofix® is a single use device.

Disposal

Disposal acc. to local guidelines and/or clinical protocols.

Storage and handling conditions

Storage conditions:

Please refer to the information on the packaging.

Notice to the user

If, during the use of this product or as a result of its use, a serious incident has occurred, please report it to the manufacturer and/or its authorised representative and to your national authority.

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B.Braun (Thailand) Ltd.,
B. Braun Melsungen AG