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EN
MANUFACTURER'S INSTRUCTIONS AND INFORMATION - INSTRUCTIONS
FOR USE OF THE PERSONAL PROTECTIVE EQUIPMENT

easyCARE latex PF Latex powder-free examination gloves
REF: RLBRRPPP07

RR - size, PPP - number of pieces per pack

Product features

- raw material: 100% natural rubber latex
- colour: white and cream
- cuff: beaded
- powder: powder-free, the content of powder residue below 2mg per glove
- AQL: 1.0
- internal surface: chlorinated
- outer surface: micro-textured on entire glove with additional texture on the fingertips
- allowed for contact with food: in accordance with Regulation (EC) 1935/2004 and Commission Regulation (EC) 10/2011
- intermediate packaging: 100 pcs with a colour-coded packaging for easy size recognition
- size: XS-XL
- single-use
- non-sterile
- shelf life: 5 years

Intended use

The device is intended for patient diagnostics, providing a barrier to the transmission of microorganisms between the operator and the patient, minimizing the risk of cross-contamination.

Product classification

- Class I Medical Device according to MDR 2017/745 (EU)

• Personal Protective Equipment Cat.III according to Regulation (EU) 2016/425
EU Type Examination and Ongoing Conformity by Notify Body:
SATRA Technology Europe Ltd, Bracetown Business Park Clonee, D15 YN2P,
Ireland (Notified Body: 2777)

Compliance with standards:

EN ISO 21420:2020+A1:2024, EN ISO 374-1:2016+A1:2018, EN ISO 374-2:2019,
EN 16523-1:2015+A1:2018, EN ISO 374-4:2019, EN ISO 374-5:2016, EN 455-
1:2020+A2:2024, EN 455-2:2024, EN 455-3:2023, EN 455-4:2009, ASTM F1671,
ASTM D3578

Contraindications

The gloves contain natural rubber latex and other ingredients used at the manufacturing stage that may cause allergic reactions, including anaphylactic shock. If these reactions occur, discontinue use of the gloves and contact a doctor. Do not use on people who are allergic to latex.

Research conducted:

Chemical permeation (EN ISO 374-1:2016+A1:2018/Type B)			Level	Medium degradation (EN ISO 374-4:2019)	
K	40%	Sodium hydroxide	6	-43,5%	Degradation levels indicate changes in puncture resistance after exposure to the triggering chemical.
M	65%	Nitric acid	1	57,0%	
P	30%	Hydrogen peroxide	1	-1,9%	
S	40%	Hydrofluoric acid	3	X	
T	37%	Formaldehyde	3	-28,3%	

The properties of the gloves when degraded by 40% hydrofluoric acid are not known, so it is important to carefully observe the physical properties of the gloves during use and to discontinue use immediately if signs of degradation occur.

X: indicates that the glove has not been tested or the test method appears to be unsuitable for the construction or material of the glove.

Levels of resistance to penetration:

Level 1: >10 min; Level 2: >30 min; Level 3: >60 min;
Level 4: >120 min; Level 5: >240 min; Level 6: >480 min

Resistance to bacteria and fungi - PASS, according to EN ISO 374-5:2016

Resistance to viruses - PASS, according to EN ISO 374-5:2016

User information:

1. Application:

These gloves are intended for single use only, on one individual during a single procedure, and should be replaced with new gloves after use.

2. Time of use:

The information on the duration of use does not reflect the actual duration of protection in the workplace and the differentiation between mixtures and pure chemicals.

3. Chemical resistance:

The chemical and penetration resistance of the glove has been assessed under laboratory conditions from samples taken from the palm only, and relates

only to the chemicals tested (this may be different if the chemical is used in a mixture).


4. Workplace protection:

- It is advisable to check that the gloves are suitable for the intended use, as workplace conditions may differ from the type test depending on temperature, abrasion and degradation.

- During use, protective gloves may provide less resistance to a hazardous chemical due to changes in physical properties. Movement, chafing, friction, degradation due to chemical contact, etc. can significantly reduce the actual time of use. For corrosive chemicals, degradation may be the most important factor to consider when selecting chemical-resistant gloves.

- Before use, check that the gloves are free from defects or damage, and that the packaging is not damaged (torn).

Putting on and taking off gloves:

Sizes:

X-Small	5-6
Small	6-7
Medium	7-8
Large	8-9
X-Large	9-10

Additional information:

- Mechanically packed in layers within the package.
- Additional information (including Declaration of Conformity) on our website www.zarys.com
- easyCARE gloves protect against hazardous substances and mixtures as well as harmful biological agents.

Special storage or handling conditions

- The product has a shelf life of 5 years from the date of manufacture.
- Store in a cool, dry, well-ventilated and clean place at: 5°C-35°C away from direct sunlight. Keep away from sources of ozone or fire. Gloves should be stored in the original packaging. Improper storage reduces the protective properties of the gloves.

Disposal

Used devices may pose a potential biohazard. They should be disposed of in accordance with local medical waste regulations.

The easyCARE gloves are classified as a Class I medical device in accordance with Regulation 2017/745 (Annex VIII) and a Category III PPE in accordance with EP and Council Regulation (EU) 2016/425 (Annex I), however, this document does not constitute instructions for use of the medical device. As a Class I medical device, the device can be safely used without instructions for use and based on the provisions of Annex I p. 23.1(d) of Regulation (EU) 2017/745, instructions for use for the medical device are not required. At the same time, these instructions in no way constitute instructions for use for a medical device and must not be used for that purpose.

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