previoMedical Set of disposable hospital bedding, non-sterile

Manufacturer: ZARYS International sp. z o.o. ul. Ziemska 44, 41-803 Zabrze, Poland

Intended use

- Disposable medical device used to cover and protect hospital beds, quilts, and pillowcases that patients come into contact with. It is used to prevent the spread of infectious agents and minimize the presence of bacterial reservoirs in the hospital environment that can be easily transmitted between patients through contact with a common surface, thereby reducing the risk of wound infection and preventing disease.
- For the protection of hospital beds, couches, trolleys in hospital wards, isolation rooms, post-operative wards, emergency departments for patients with a recognised infectious disease or infection and in the case of a suspected outbreak of an alert pathogen in biological material. According to current recommendations for preventing the spread of nosocomial infections and minimising the presence of bacterial reservoirs in the hospital environment, a set of disposable medical bedding is recommended to be used routinely in the event of patient isolation. This standard of practice protects the patient from cross-infection, which originates from the transmission to the patient of an infectious agent, bacteria, viruses and/or fungi from another patient, which has a direct impact on the health and life of patients, particularly vulnerable patients, critically ill patients, those with extensive wounds as well as those with low immune resistance.
- In order to eliminate the undesirable phenomenon of the release of "dust" particles occurring during the use of reusable bedding sets made of cotton, a polypropylene/SMS non-dusting non-woven fabric was used in the production of the product. These non-woven fabrics minimise the possibility of dust and "dust" penetrating into the wound, significantly reducing the risk of infection posing a threat to the patient's health and life.
- It also reduces the transport of contaminated reusable bedding sets to the laundry and the implementation of additional epidemiological procedures, as well as minimising the risk of indirect persons coming into contact with a potential pathogen.

Ethylene oxide sterilisation is required for use in aseptic conditions.

Indications for use

- The need to have procedures in place to ensure that patients are protected from infections and infectious diseases.

Place of application

- Hospital wards, isolation rooms, post-operative wards, emergency rooms, emergency departments

Target group

- Patients irrespective of gender and age, isolated patients, patients particularly susceptible to infection, seriously ill patients, patients with extensive wounds as well as those with low immune resistance, patients undergoing surgery

Intended user

- Professional user - medical personnel

Clinical benefits

- Reducing the transmission of contaminants and/or potential infectious agents that may pose a health risk to the patient
- Protecting the patient from infection by an alarm pathogen during medical procedures in hospital wards
- Minimising pathogen transmission during patient isolation
- Reducing the risk of infection for auxiliary staff involved in the transport and laundering of contaminated linen

Maximum duration of use

- The use time of the product should not exceed 24 hours

Preparation of the product for use

Disposable, non-sterile device. If used under aseptic conditions, the device must be sterilised with ethylene oxide.

Sterilisation with ethylene oxide in a process validated in accordance with the current standard EN ISO 11135 is recommended. The sterilisation process should take place in sterilisers meeting the requirements of the current standard EN 1422.

The sterilisation process requires repackaging of a single unit in accordance with the hospital procedure in packaging conforming to the applicable standards EN ISO 11607-1 and EN ISO 11607-2. Microbiological cleanliness level of the product max. 300 CFU/100 cm².

STERILISATION CONDITIONS:

- pre-conditioning: 99 min
- conditioning: 60 min
- air removal: 11 min
- pressure time: 10 min
- EO injection: 20 min
- EO removal: 33 min
- rinsing: 179 min

STERILISATION PROCESS PARAMETERS:

- temperature (°C): 47-57
- pressure (KPa): -15 ±2
- humidity (%RH): 30-70
- exposure time (min): 300

After sterilisation, the device must be quarantined according to the hospital procedure for ethylene oxide sterilisation.

Application instructions

1. It is necessary to select the size of the individual pieces on which the bedding will be placed.
2. Check the packaging for damage.
3. Open the packaging, then remove the individual components hygienically and check for damage before use.
4. Place the unfolded sheet over the mattress, pillowcase and pillowcase, etc. Product fitted with tabs to secure and hold properly.
5. If any part of the set is damaged, contaminated, it must be replaced with a new one.
6. Disposable device, replace after each patient.

Product notes and warnings

- Use the product in accordance with its intended use
- Use immediately after taking out of the packaging
- Hygiene rules must be observed when using the product
- Do not use the product if its packaging is damaged
- Do not use the product if it is damaged or contaminated
- Single-use device, cannot be reused, reuse endangers the health and life of the patient and may lead to an increased risk of infection or cross-contamination
- Do not use the product beyond the safe use date stated on the packaging
- If symptoms of irritation or allergic reaction occur, remove the device and discontinue its use
- Storage in the original packaging reduces the risk of microbiological contamination

Special storage or handling conditions

- The product has a shelf life of 5 years from the date of manufacture
- Store in accordance with the conditions specified on the product packaging: Keep away from sunlight, protect from moisture

Disposal

Used devices may pose a potential biohazard. Devices should be disposed of in accordance with local medical waste regulations.

Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State.

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