EN
INSTRUCTIONS FOR USE
previoMedical Procedure gown, with knitted cuffs, non-woven, non-sterile
previoMedical Procedure gown, with elastic cuffs, non-woven, non-sterile

Manufacturer: ZARYS International sp. z o.o. ul. Ziemska 44, 41-803 Zabrze, Poland

Intended use

- Disposable medical device intended for use by medical personnel in hospital wards and emergency rooms during medical procedures that do not require sterile conditions. It is used as a barrier to minimize the transfer of dust particles and exfoliated skin, which are carriers of pathogens, from the surface of the staff's body to the patient, thus reducing the risk of wound infection and preventing disease. It is also used during procedures on patients with low immune resistance, who are particularly vulnerable to infection, and on patients with extensive wounds, due to the opening of a gateway for potential pathogen invasion.
- Indicated for use during aseptic procedures such as dressing changes, bladder catheterisation, suctioning of respiratory secretions, blood sampling for tests, administration of medication.
- Used to protect the patient from cross-infection, which originates from the transmission to the patient of an infectious agent, bacteria, viruses and/or fungi from another patient.
- Recommended during diagnostic examinations, e.g. gastroscopy, colonoscopy, bronchoscopy
- In order to eliminate the undesirable phenomenon of the release of "dust" particles occurring during the use of reusable aprons made of cotton, a non-dusting polypropylene non-woven fabric is used in the manufacture of the product. This non-woven fabric minimises the possibility of dust and "dust" entering the wound, significantly reducing the risk of infection posing a threat to the patient's health and life.

Ethylene oxide sterilisation required for use in aseptic conditions

Indications for use

- The need to protect the patient due to the opening of the floodgates of potential invasion by a pathogenic agent

Place of application

- Isolation rooms, recovery rooms, treatment wards, outpatient clinics, EDs, treatment rooms, intensive care units.

Target group

- Patients irrespective of gender and age, patients with extensive wounds, patients with low immunological immunity, patients requiring invasive diagnostic tests

Intended user

- Professional user - medical personnel

Clinical benefits

- Protecting the patient from infection with extensive wounds.
- Protecting the patient from infection with low immune resistance
- Protecting the patient from cross-infection during invasive diagnostic examinations such as gastroscopy, colonoscopy, bronchoscopy

due to the opening of the floodgates of a potential invasion by a pathogenic agent

Maximum duration of use

- The use time of the product should not exceed 24 hours

Preparation of the product for use

Disposable, non-sterile device. If used under aseptic conditions, the device must be sterilised with ethylene oxide.

Sterilisation with ethylene oxide in a process validated in accordance with the current standard EN ISO 11135 is recommended. The sterilisation process should take place in sterilisers meeting the requirements of the current standard EN 1422.

The sterilisation process requires repackaging of a single unit in accordance with the hospital procedure in packaging conforming to the applicable standards EN ISO 11607-1 and EN ISO 11607-2.

STERILISATION CONDITIONS:

- pre-conditioning: 99 min
- conditioning: 60 min
- air removal: 11 min
- pressure time: 10 min
- EO injection: 20 min
- EO removal: 33 min
- rinsing: 179 min

STERILISATION PROCESS PARAMETERS:

- temperature (°C): 47-57
- pressure (KPa): -15 ±2
- humidity (%RH): 30-70
- exposure time (min): 300

After sterilisation, the device must be quarantined according to the hospital procedure for ethylene oxide sterilisation.

Application instructions

1. The size of the apron should be matched to the dimensions of the wearer.
2. Check the packaging for damage.
3. Open the packaging, then remove the apron under hygienic conditions and check for damage before use.
4. Wear the apron so that the front of the body is completely covered and the slit is on the user's back side.
5. Tie the apron around the waist and around the neck - use an extra person if necessary
6. If damaged, contaminated, the apron should be replaced with a new one.
7. Disposable device, replace after contact with each patient.

Product notes and warnings

- Use the product in accordance with its intended use
- Use immediately after taking out of the packaging
- Hygiene rules must be observed when using the product
- Do not use the product if its packaging is damaged.
- Do not use the product if it is damaged or contaminated
- Single-use device, cannot be reused; reuse endangers the health and life of the patient and may lead to an increased risk of infection or cross-contamination
- Do not use the product beyond the safe use date stated on the packaging
- If symptoms of irritation or allergic reaction occur, remove the device and discontinue use
- Storage in the original packaging reduces the risk of microbiological contamination

Special storage or handling conditions

- The product has a shelf life of 5 years from the date of manufacture
- Store in accordance with the conditions specified on the product packaging: Keep away from sunlight, protect from moisture

Disposal

Used devices may pose a potential biohazard. Devices should be disposed of in accordance with local medical waste regulations. Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State.

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